

Usono ProbeFix System -

For use with all brands of cardiac ultrasound transducers

Echocardiography is one of the essential imaging techniques for the diagnosis of cardiovascular diseases.

Many specialists face significant difficulties when acquiring a stable image, scanning obese patients, and working under pressure.


Imaging during exercise is particularly challenging due to the patient's vigorous chest movements. Furthermore, specialists increasingly experience pain while scanning. Redundant arm abduction and high muscle activity cause severe soreness in shoulders, arms and wrists.

In collaboration with multiple specialists, an accessory to attach the probe to the chest has been developed. The ProbeFix is a modular tool that attaches an existing cardiac ultrasound probe to the body.

The probe can still be physically tilted and rotated at any time. With 3D/4D probes with digital rotation the workflow is even easier and more hands-off.

The stable and lengthy fixation provides continuous and better reproducible images.

The window stays the same, even during exercise on an ergometer or a treadmill. Ergonomics are improved thanks to the hands-off workflow, which can eventually prevent repetitive strain injuries.



Features:

Hygienic Straps

- Fully disposable high-quality straps ensure optimal hygiene.

Modular Probe Adapter

- Different Probe Adapters guarantee a perfect fit of each ultrasound probe.

Silicon Feet

- Silicon feet combine perfect grip on the body with great comfort for patients.

Sphere Mechanism

- Carefully designed mechanism facilitates effortless movement into the right position.

Perfect Stable Position

- A patented mechanism ensures a perfect stable position of the probe, for as long as required.

Easy To Clean

- Take apart easily to clean all different components.

Usono

Usono ProbeFix System

Ref. no.	Qty/box	Description
P-1001	1	ProbeFix base assembly
P-2011	1	ProbeFix adaptor for Phillips S5-1
P-2012	1	ProbeFix adaptor for Phillips S4-2
P-2013	1	ProbeFix adaptor for Phillips X5-1
P-2021	1	ProbeFix adaptor for GE M5Sc
P-2022	1	ProbeFix adaptor for GE M3Sc
P-2031	1	ProbeFix adaptor for Siemens P4-2
P-2041	1	ProbeFix adaptor for Esaote PA-230E
P-9001	10	ProbeFix straps (small)
P-9002	20	ProbeFix straps (medium)
P-9003	10	ProbeFix straps (large)
P-9004	10	ProbeFix straps (extra large)
P-9010	20	ProbeFix straps (shoulder)

DS-10-20 Version 3

Contact/Place order: Mermaid Medical A/S · Frydensbergvej 25 · DK-3660 Stenløse · www.mermaidmedical.com · customer@mermaidmedical.com



Denmark
P +45 47 10 85 70
F +45 47 10 85 71

Sweden
P +46 31 755 41 28
F +46 31 755 45 99

Finland
P +358 923 113238
F +358 923 113239

Norway/Iceland
P +47 21 53 88 94
F +47 21 53 88 95

United Kingdom
P +44(0)1704 560493
F +44(0)1704 551182

Spain
P +34 9253 33159
F +34 9253 33159
customer@mermaidmedical.es

Netherlands
P +31 314 63 2070
F +31 314 63 1685
info@mermaidmedical.nl

Germany
P +49 40 8518665-0
F +49 40 8518665-1
info@mermaidmedical.de