



## Scholten Novatome™

A biopsy design proven by over 30 years of successful use. The same sharp and steerable advantages in a new disposable format.

Scholten Novatome™ is specially designed to acquire maximum specimen size for more accurate pathology results.

The two-ringed handle easily adapt to any thumb position. The handle provides a solid instrument feel. With a reputation for its ease of positioning and heart rhythm feedback, Scholten Bioptomes have earned high praise.

The Novatome™ Single use, disposable unit eliminates the risk of contamination. Scholten model is the best choice for safety and precision. Clearly, with a pedigree of more than 30 years of clinical experience, the new disposable model Novatome™ will continue to set the standard worldwide.



Biopsy Forceps Models				
Part number	French size	Length	Curved tip	Jaw cup dia.
SU104-50	6	50 cm	Yes	2.0 mm
SU104-100	6	100 cm	Yes	2.0 mm
SU103-50	7	50 cm	Yes	2.3 mm
SU103-100	7	100 cm	No	2.3 mm
SU103-100 PC	7	100 cm	Yes	2.3 mm
SU102-50	8	50 cm	Yes	2.7 mm
SU101-50	9	50 cm	Yes	3.0 mm

- Features:**
- Perfect balance of flexibility and stiffness offering precise navigation and positioning
  - Provides precise cutting edge alignment for every cut
  - Formable distal curve
  - Single action cutting jaw (Single moving jaw) provides superior cutting force
  - Familiar locking finger ring actuation handle
  - Unique hollow ground cutting edge for superior sharpness
  - Super precision machining

DS-4-04 Version 3

Contact/Place order: Mermaid Medical A/S · Frydensbergvej 25 · DK-3660 Stenløse · www.mermaidmedical.com · customer@mermaidmedical.com



**Denmark**  
P +45 47 10 85 70  
F +45 47 10 85 71

**Sweden**  
P +46 31 755 41 28  
F +46 31 755 45 99

**Finland**  
P +358 923 113238  
F +358 923 113239

**Norway/Iceland**  
P +47 21 53 88 94  
F +47 21 53 88 95

**United Kingdom**  
P +44(0)1704 560493  
F +44(0)1704 551182

**Spain**  
P +34 9253 33159  
F +34 9253 33159  
customer@mermaidmedical.es

**Netherlands**  
P +31 314 63 2070  
F +31 314 63 1685  
info@mermaidmedical.nl

**Germany**  
P +49 40 8518665-0  
F +49 40 8518665-1  
info@mermaidmedical.de